

HCP-7000
CHART PROJECTOR

HUVITZ
CHART PROJECTOR





HCP-7000 CHART PROJECTOR

- Semi-permanent lifetime LED light source
- Brighter and clearer than conventional bulb light
- Attractive Streamlined Shape and Blue Indicator
- RC-R Type(Russian Chart) added

SPECIFICATIONS

Chart	41 charts 34 masks, Red / Green & Polarization Filters
Distance of projection	2.5-8m(Optional: 8m screen)
Chart Rotation Speed	Average 0.15 sec
Chart magnification	30X (When 5m distance)
Power Saving	Automatic Lamp Off (10 min)
Program	2 programs with a maximum of 30 charts each
Tilt Angle	15°
Power Supply	100-120V@50Hz : 1.0A 200-240V@60 Hz : 0.5A
Lamp	LED 4W
Fuse	250V, T3.15AL
Dimension	270(W) x 182(D) x 230(H)mm
Weight	3.44kg (Body: 3.06kg + Stand: 0.38kg)

Designs and details can be changed without prior notice for the purposes of improvement.